



Food and Drug Branch – Device Recalls

CALIFORNIA DEVICE RECALL INFORMATION SHEET

Resmed Ltd. Recalls Astral 100 And Astral 150 Ventilators For Sudden Power Loss

Recall Date	Product Description	Recalling Firm	Recall Reason
11/1/2023	Astral 100 And Astral 150 Ventilators	RESMED LTD. Bella Vista, Australia	If ventilator is on internal battery, not intended to serve as a primary power source, low/critically low battery alarms will sound, but a fault leads to sudden power loss. If power fails, then ventilation stops and a Total Power Failure (TPF) alarm should sound, but it's powered by a supercapacitor, which degrades over time, which may cause TPF alarm to sound for less than 2 minutes or not at all

Recall Class	Product Identification	Distribution	Affected Dates
I	Astral 100 And Astral 150 Ventilators All devices manufactured between 2013 and 2019 Astral 100, UDI: 00619498003259 Astral 150, UDI: 00619498003266	16634 Units Nationwide including California	13 September 2023 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

