Screening Information System Project Genetic Disease Screening Program – Prenatal Screening

View / Refresh Interps

Screen Negative

Screen Positive *or*Further Information Needed

<u>Screen Negative</u> results are immediately available once you have saved NT data in SIS.

Patients with a first trimester preliminary Screen Negative result should submit a second trimester blood sample to obtain Full Integrated Screening and a refined risk assessment. Initial Screen Positive results require validation of all interpretation factors before a result can be finalized

Missing interpretation
factors will require validation
before an interpretation can
be created

In both cases, you will see a message asking you to contact the Coordinator. All interpretation factors must be verified before the interpretation can be recalculated and finalized.

First trimester interpretation factors include: blood collection date, race/ethnicity, patient weight, ovum donor age (if relevant), diabetic status, and smoking status.

Screening Information System Project Genetic Disease Screening Program – Prenatal Screening



You wish to provide your patient with results

You do not wish to proceed with interpretation verification

- 1. Call the Coordinator in order to verify the interpretation factors for your patient.
 - 2. Wait as the Coordinator reinterprets the case in SIS. Click on the "View/Refresh Interps" button to see the recalculated results.
- 3. Click on the "Print Interps" button to provide the patient with a copy of the results, as well as information about her follow-up options authorized through the PNS Program.

- 1. You can exit from SIS. The Coordinator will be alerted through SIS that the interpretation factors must be reviewed, and will contact the referring clinician.
- 2. The referring clinician will receive the final interpretation result and be responsible for communicating the results to the patient.