



State of California—Health and Human Services Agency  
California Department of Public Health



Food and Drug Branch – Device Recalls

**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**Abbott Diabetes Care Recall Components of FreeStyle Libre Flash Glucose Monitoring Systems for Lithium Battery Problem**

Recall Date	Product Description	Recalling Firm	Recall Reason
4/12/2023	<b>FreeStyle Libre Reader</b> REF: 71525-01, 71701-01, 71739-01 a component of the FreeStyle Libre Flash Glucose Monitoring System	<b>Abbott Diabetes Care, Inc.</b> Alameda, California	Lithium-ion batteries in glucose monitor- ing system readers may swell, overheat, or pose fire hazard. Only use the charg- ing cable/power ada- pter included with re- aders. Do not use reader if it is damag- ed, cracked, swell- ing, too hot to hold, it will no longer hold charge, or it fails Re- ader Test. Do not use other chargers. User Manual update includes instructions on reader safe stor- age, charging, use.
4/12/2023	<b>FreeStyle Libre Reader</b> REF: 71936-01, 71937-01, 71938-01 a component of the FreeStyle Libre 14 Day Flash Glucose Monitoring System	<b>Abbott Diabetes Care, Inc.</b> Alameda, California	Lithium-ion batteries
4/12/2023	<b>FreeStyle Libre Reader</b> REF: 71951-01, 71952-01, 71953-01 a component of the FreeStyle Libre 2 Flash Glucose Monitoring System	<b>Abbott Diabetes Care, Inc.</b> Alameda, California	Lithium-ion batteries

Recall Class	Product Identification	Distribution	Affected Dates
I	<b>FreeStyle Libre Flash Glucose Monitoring System Reader</b> REF/UDI-DI/Lot: 71525-01/357599000219/02Y702D -02Y304D ; 71701-01/357599000226/01Y702D -02Y317D; 71739-01/357599000240/03Y754D -01Y671D	336085 units Nationwide	February 2023 and prior
I	<b>FreeStyle Libre 14 Day Flash Glucose Monitoring System Reader</b> REF/UDI-DI/Lot: 71936-01/357599000981/01C083D -03Z057D; 71937-01/357599000974/02B099D -02Y763D; 71938-01/357599002008/04C216D -03Y760D	1373311 units Nationwide	February 2023 and prior
I	<b>FreeStyle Libre 2 Flash Glucose Monitoring System Reader</b> REF/UDI-DI/Lot: 71951-01/357599805005/02C313D -02A279D; 71952-01/357599804008/01B179D -04A205D; 71953-01/357599803001/06B074F -03A198D	1534532 units Nationwide	February 2023 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch  
 MS 7602 • P.O. Box 997435 • Sacramento, CA 95899-7435  
 (916) 650-6500 • (916) 650-6650 FAX  
 Internet Address: [www.cdph.ca.gov](http://www.cdph.ca.gov)

