



**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**Invacare Corporation Recalls Dolomite Gloss Rollator For Premature Failure Of The Seat During Use, The Plastic Eyelets Of The Folding Mechanism Of The Device's Seat Can Break**

Recall Date	Product Description	Recalling Firm	Recall Reason
8/30/2023	<b>Invacare Corporation</b> DOLOMITE GLOSS 680 ROLLATOR	<b>Invacare Corporation</b> Elyria, Ohio	Premature failure of the seat during use
8/30/2023	<b>Invacare Corporation</b> DOLOMITE GLOSS 600 ROLLATOR	<b>Invacare Corporation</b> Elyria, Ohio	Premature failure of the seat during use
8/30/2023	<b>Invacare Corporation</b> DOLOMITE GLOSS 520 ROLLATOR	<b>Invacare Corporation</b> Elyria, Ohio	Premature failure of the seat during use
8/30/2023	<b>Invacare Corporation</b> DOLOMITE GLOSS 450 ROLLATOR	<b>Invacare Corporation</b> Elyria, Ohio	Premature failure of the seat during use

Recall Class	Product Identification	Distribution	Affected Dates
II	<b>Invacare Corporation</b> DOLOMITE GLOSS 680 ROLLATOR Model Number(s): 1654807, Serial number 20BHL0002 to 21IHL1910	586 devices Nationwide	June 2023 and prior
II	<b>Invacare Corporation</b> DOLOMITE GLOSS 600 ROLLATOR Model Number(s): 1655081, Serial number 20BHL0002 to 21IHL1910	586 devices Nationwide	June 2023 and prior

II	<b>Invacare Corporation</b> DOLOMITE GLOSS 520 ROLLATOR Model Number(s): 1655082, Serial number 20BHL0002 to 21IHL1910	586 devices Nationwide	June 2023 and prior
II	<b>Invacare Corporation</b> DOLOMITE GLOSS 450 ROLLATOR Model Number(s): 1655083, Serial number 20BHL0002 to 21IHL1910	586 devices Nationwide	June 2023 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

