



Food and Drug Branch – Device Recalls

CALIFORNIA DEVICE RECALL INFORMATION SHEET

Olympus Corporation Recall Evis Bronchovideoscope Olympus Bf Type For Risk Of Endobronchial Combustion

Recall Date	Product Description	Recalling Firm	Recall Reason
11/15/2023	Bronchovideoscope Olympus Bf Type Model Numbers BF-1T150, BF-1T60, BF-P150 BF-P60 and EVIS EXERA BRONCHOVIDEOSCOPE OLYMPUS BF TYPE Model Number XT160	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen [and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.
11/15/2023	Evis Exera li Bronchovideoscope Olympus Bf Type Model Number Q180-AC	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen [and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.
11/15/2023	EVIS EXERA II BRONCHOVIDEOSCOPE OLYMPUS BF TYPE Model Numbers BF-1T180, BF-1TQ180, BF-P180	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen

			[and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.
11/15/2023	EVIS EXERA III BRONCHOVIDEOSCOPE Model Numbers BF-1TH190, BF-H190, BF-Q190 and BRONCHOVIDEOSCOPE OLYMPUS BF-Q170 Model Number BF-Q170	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen [and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.
11/15/2023	EVIS EXERA b BRONCHOVIDEOSCOPE OLYMPUS Model Number BF-XT190	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen [and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.
11/15/2023	EVIS EXERA III BRONCHOVIDEOSCOPE OLYMPUS Model Number BF-P190	OLYMPUS CORPORATION OF THE AMERICAS Center Valley, Pennsylvania	There is a risk of endobronchial combustion if high-frequency cauterization is performed while supplying oxygen [and/or] the electrode section of the electrosurgical accessory is too close to the distal end of the endoscope.

Recall Class	Product Identification	Distribution	Affected Dates
I	Bronchovideoscope Olympus Bf Type Model BF-1T150 UDI-DI: 4953170308185; Model BF-1T60 UDI-DI: 4953170339264; Model BF-P150 UDI-DI: N/A; Model BF-P60 UDI-DI: 4953170339196; Model BF-XT160 UDI-DI: N/A	15942 Units Nationwide including California	12 October 2023 and prior
I	EVIS EXERA II BRONCHOVIDEOSCOPE OLYMPUS BF TYPE Model BF-Q180-AC UDI-DI: N/A	1154 Units Nationwide including California	12 October 2023 and prior
I	EVIS EXERA II BRONCHOVIDEOSCOPE OLYMPUS BF TYPE Model BF-1T180 UDI-DI: N/A; Model BF-1TQ180 UDI-DI: N/A; Model BF-P180 UDI-DI: N/A;	13950 Units Nationwide including California	12 October 2023 and prior
I	EVIS EXERA III BRONCHOVIDEOSCOPE Model BF-1TH190 UDI-DI: 4953170335181; Model BF-H190 UDI-DI: 4953170335174; Model BF-Q170 UDI-DI: 4953170342912 Model BF-Q190 UDI-DI: 4953170335198	32211 Units Nationwide including California	12 October 2023 and prior
I	EVIS EXERA b BRONCHOVIDEOSCOPE OLYMPUS Model BF-XT190 UDI-DI: 4953170402470	719 Units Nationwide including California	12 October 2023 and prior

I	EVIS EXERA III BRONCHOVIDEOSCOPE OLYMPUS Model BF-P190 UDI-DI: 4953170342110	3948 Units Nationwide including California	12 October 2023 and prior
---	---	--	------------------------------

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch
MS 7602 ● P.O. Box 997435 ● Sacramento, CA 95899-7435
(916) 650-6500 ● (916) 650-6650 FAX
Internet Address: www.cdph.ca.gov

