



California Medical Device Recall Information



Recall Name

Abbott Diabetes Care Recalls FreeStyle Insulinx Blood Glucose Meters Due to Incorrect Test Results

Recall Date	Product Description	Recalling Firm	Recall Reason
04/15/13	FreeStyle Insulinx Blood Glucose Meters are used to measure, log, and generate reports of insulin and blood glucose test results.	Abbott Diabetes Care Alameda, CA	<i>At extremely high blood glucose levels (1024 mg/dL and above), the meter will report an incorrect test result.</i>
Recall Class	Product Identification	Distribution	Affected Dates
TBD	All units manufactured from product launch through April 1, 2013.	CA , nationwide	Distributed between April 18, 2012 and April 1, 2013.

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/Safety/Recalls/ucm348391.htm?source=govdelivery>