



# California Medical Device Recall Information



## Recall Name

### Alere Recalls INRatio® and INRatio®2 PT/INR Monitor System Due to Incorrect Results

Recall Date	Product Description	Recalling Firm	Recall Reason
Updated: 12/08/14  Initial: 04/16/14	<ul style="list-style-type: none"><li>• INRatio®2 PT/INR Monitors</li><li>• INRatio® PT/INR Monitors</li><li>• INRatio® PT/INR Test Strips</li></ul>	<b>Alere, Inc.</b> San Diego, CA	<i>In certain cases an INRatio® and INRatio®2 PT/INR Monitor System may provide inaccurate INR results.</i>  <i>Issues can arise if the patient has certain medical conditions or if the System instructions in the labeling for performing the test are not followed.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	<b>All INRatio® PT/INR test strips and monitors.</b>  <a href="#">Product Photos</a>	<b>CA</b> , nationwide	Test strips distributed from:  September 4, 2013 to December 4, 2014.  Monitors were distributed from:  April 2008 to December 3, 2014.

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm426166.htm>