



California Medical Device Recall Information



Recall Name

Baxter Recalls VASCU-Guard Peripheral Vascular Patch Due to Potential for Incorrect Patch Orientation

Recall Date	Product Description	Recalling Firm	Recall Reason
05/02/15	VASCU-Guard Peripheral Vascular Patch	Baxter International, Inc. Deerfield, IL	<i>Baxter received customer complaints of difficulty in distinguishing the smooth from rough surface of the VASCU-GUARD patch.</i> <i>Incorrect orientation of the patch with the rough side toward the bloodstream may increase the risk of vessel thrombosis and/or embolism.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	<ul style="list-style-type: none">Vascu-Guard TS, 1x6cm Product code 1504026Vascu-Guard TS, 0.8x8cm Product code 1504028Vascu-Guard TS, 1x10cm Product code 1504030Vascu-Guard TS, 2x9cm Product code 1504032	CA , nationwide	Manufacturing dates: January 19, 2015 to May 1, 2015 Distribution dates: March 16, 2015 to May 1, 2015

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm449825.htm>