



California Medical Device Recall Information



Recall Name

Covidien Recalls Monoject Prefill Flush Syringes Due to Sterility Assurance and Quality Control Concerns

Recall Date	Product Description	Recalling Firm	Recall Reason																														
8/16/13	<p>Monoject prefill flush 12mL Syringes (multiple fill volumes affected):</p> <ul style="list-style-type: none">• Monoject 0.9% Sodium Chloride Flush Syringe• Monoject 10 Units/mL Heparin Lock Flush• Monoject 100 Units/mL Heparin Lock Flush	Covidien Mansfield, MA	<i>Possible risk that syringes were not subjected to autoclave sterilization.</i> <i>Some syringes also have mismatched caps, labels, and wrapper.</i>																														
Recall Class	Product Identification	Distribution	Affected Dates																														
I	<table><tr><th><u>Product ID</u></th><th><u>Lot #*</u></th></tr><tr><td>8881570121</td><td>13A0084N</td></tr><tr><td>“</td><td>13A0094</td></tr><tr><td>“</td><td>13B0364</td></tr><tr><td>“</td><td>13C0504</td></tr><tr><td>“</td><td>13C0514</td></tr><tr><td>8881570123</td><td>13A0084N</td></tr><tr><td>8881570125</td><td>13A0084N</td></tr><tr><td>8881580121</td><td>13A0084N</td></tr><tr><td>8881580123</td><td>13A0084N</td></tr><tr><td>8881580125</td><td>13A0084N</td></tr><tr><td>8881590121</td><td>13A0084N</td></tr><tr><td>8881590123</td><td>13A0084N</td></tr><tr><td>8881590125</td><td>13A0084N</td></tr><tr><td>“</td><td>13D0824N</td></tr></table> <p>* Lot numbers can be found on the shipper case, carton and individual syringes.</p>	<u>Product ID</u>	<u>Lot #*</u>	8881570121	13A0084N	“	13A0094	“	13B0364	“	13C0504	“	13C0514	8881570123	13A0084N	8881570125	13A0084N	8881580121	13A0084N	8881580123	13A0084N	8881580125	13A0084N	8881590121	13A0084N	8881590123	13A0084N	8881590125	13A0084N	“	13D0824N	CA , nationwide	Covidien alerted customers of the recall by letter on 8/16/13.
<u>Product ID</u>	<u>Lot #*</u>																																
8881570121	13A0084N																																
“	13A0094																																
“	13B0364																																
“	13C0504																																
“	13C0514																																
8881570123	13A0084N																																
8881570125	13A0084N																																
8881580121	13A0084N																																
8881580123	13A0084N																																
8881580125	13A0084N																																
8881590121	13A0084N																																
8881590123	13A0084N																																
8881590125	13A0084N																																
“	13D0824N																																

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/Safety/Recalls/ucm365577.htm>