



California Medical Device Recall Information



Recall Name

**Dräger Recalls the PS500 Battery Power Supply
Used With Evita V500 and Babylog VN500 Ventilators
Due to Software Issue Causing Short Battery Run Time**

Recall Date	Product Description	Recalling Firm	Recall Reason
<p>Initial: 7/13/15</p> <p>Expanded: 12/01/15</p>	<p>Option PS500 Battery Power Supply</p> <p>Used with:</p> <ul style="list-style-type: none"> • Evita V500 Ventilator • Babylog VN500 Ventilator 	<p>Dräger Medical, Inc. Telford, PA</p>	<p><i>A software issue causes shorter than expected battery run times.</i></p> <p><i>The software issue with the power supply also prevents the appropriate alarm from sounding 5 minutes before the battery runs out of power and the device shuts down.</i></p>
Recall Class	Product Identification	Distribution	Affected Dates
<p>I</p>	<p>All Units Affected</p> <p>Catalog Numbers: 8416400 - 8417400</p>	<p>CA, nationwide</p>	<p>Manufactured and Distributed between:</p> <p>June 1, 2011 and October 30, 2015</p>

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm480135.htm>