



California Medical Device Recall Information



Recall Name

**GE Healthcare Recalls Giraffe Incubator with Servo Oxygen Control System
and Giraffe OmniBed
Due to Inaccurate Temperature and Oxygen Settings**

Recall Date	Product Description	Recalling Firm	Recall Reason
01/08/13	Giraffe Incubator with Servo Control Oxygen System and Giraffe OmniBed.	GE Healthcare, LLC. Waukesha, WI	<i>Manually entered oxygen and temperature settings may be overwritten during power on or recovery from a power failure.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	Units manufactured from 08/01/12 to 12/20/12. Giraffe Incubators without a Servo Control Oxygen System are not impacted.	CA , nationwide	Manufactured from 08/01/12 to 12/20/12.

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/Safety/Recalls/ucm349044.htm?source=govdelivery>