



California Medical Device Recall Information



Recall Name

Hospira Recalls Intravascular Administration Sets Due to a Potential of Puncture

Recall Date	Product Description	Recalling Firm	Recall Reason
Initial: 4/01/13 Updated: 8/13/13	Intravascular Administration Sets	Hospira, Inc. Lake Forest, IL	<i>Potential for the piercing pin on the Hospira Blood Sets to puncture the outer wall of non-ISO-compliant blood bags.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	Hospira Blood Sets Listed Products Recalled: <ul style="list-style-type: none">• 14200-04-28 Secondary Blood Set;• 14203-04-28 Blood Set;• 14206-04-28 Y-Type Blood Set;• 14207-04-28 Blood Set;• 14210-04-28 Plum Blood Set;• 14211-04-28 Plum Blood Set;• 14212-04-28 Plum Y-Type Blood Set;• 14217-04-28 Y-Type Blood Set• 14219-04-28 Y-Type Blood Set	CA , nationwide	Distributed from: July 2011 through February 2013

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm364988.htm>