



California Medical Device Recall Information



Recall Name

Hospira Recalls GemStar™ Docking Stations Due to Potential for Malfunction

Recall Date	Product Description	Recalling Firm	Recall Reason
05/02/14	GemStar™ Docking Station	Hospira, Inc. Lake Forest, IL	<i>The docking station may malfunction when used with a GemStar Phase 3 or Phase 4 infusion pump.</i> <i>The pump may fail to power up, or an error code and alarm may be triggered, which could stop infusion.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	List Number 13075	CA , nationwide	Hospira is retiring the GemStar™ family of infusion devices. As of July 31, 2015 , Hospira will no longer support the device.

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm395828.htm>