



California Medical Device Recall Information



Recall Name

Lumenis Limited Recalls VersaCut Morcellator Due to Potential for Air Embolisms

Recall Date	Product Description	Recalling Firm	Recall Reason
01/02/13	VersaCut Morcellators <ul style="list-style-type: none">A medical device used to break up and remove dissected tissue during surgical procedures.	Lumenis Limited Yokneam, Israel	<i>Inaccurate device labeling leading to incorrect hook up of aspirating tubing.</i> <i>Incorrect assembly could cause the introduction of air embolisms resulting in serious adverse health consequences including death.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	VersaCut Morcellators Model Numbers: <ul style="list-style-type: none">0637-245-01 (starter kit)0636-470-01 (control box)	CA , nationwide	All units manufactured since May 1998

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm340943.htm>