



California Medical Device Recall Information



Recall Name

Nidek Medical Products, Inc. Recalls Mark 5 Nuvo, Nuvo 8, and Nuvo Lite Oxygen Concentrators Due To Possible Capacitor Failure

Recall Date	Product Description	Recalling Firm	Recall Reason
5/01/12	Mark 5 Nuvo, Nuvo 8, and Nuvo Lite Oxygen Concentrators	Nidek Medical Products, Inc. Birmingham, AL	<i>Capacitor failure could result in fire hazard</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	Mark 5 Nuvo, Nuvo 8, and Nuvo Lite Oxygen Concentrators Serial Numbers affected: <ul style="list-style-type: none"> 042-10000 thru 102-09335 Serial Number is located on the rear of the concentrator	CA , nationwide	Shipped January 1, 2004 through May, 2010

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm309420.htm>