



# California Medical Device Recall Information



## Recall Name

**Stryker Orthopaedics Recalls ShapeMatch Cutting Guides  
Due to the Potential that Guides Were Not Manufactured Per Surgeon Parameters**

Recall Date	Product Description	Recalling Firm	Recall Reason
4/10/13	<b>ShapeMatch Cutting Guides</b> <ul style="list-style-type: none"><li>Custom, patient-specific surgical instrumentation for total knee replacement surgery.</li></ul>	<b>Stryker Orthopaedics / Stryker Howmedica Osteonics Corp.</b> Mahwah, NJ	<i>Defects in production software resulted in adulterated and unapproved surgical guides that when used as intended potentially results in serious adverse health consequences.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	Catalog Numbers: <ul style="list-style-type: none"><li><b>TR7100-R</b></li><li><b>TR7100-L</b></li></ul> <p>[NOTE: The <i>Triathlon Knee System</i> and <i>Triathlon</i> standard instrumentation are not affected by the recall.]</p>	<b>CA</b> , nationwide	Manufactured and distributed from:  May 2011 to November 2012.

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/Safety/Recalls/ucm347552.htm>

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm348536.htm>