



California Medical Device Recall Information



Recall Name

**Thoratec Corporation Issues Urgent Medical Device Correction Letter
for HeartMate II LVAS to Prevent Advisory Alarms and Potential Injury**

Recall Date	Product Description	Recalling Firm	Recall Reason
09/14/15	HeartMate II LVAS "Pocket" System Controller	Thoratec Corporation Pleasanton, CA	<i>Certain advisory alarms, associated with expiration of the system backup battery, have resulted in patients deciding to attempt a System Controller exchange. Injuries and deaths were reported in patients who were unable to connect to their back up System Controller in a timely manner.</i> <i>For all advisory ("yellow wrench") alarms, patients should first call their hospital contact for instructions.</i> <i>Do not attempt to replace the System Controller unless instructed by the hospital.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	All Serial Numbers that start with " PC " (i.e. "Pocket Controller")	CA , nationwide	N/A

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm462627.htm>