

Lab Data Entry Tool (LDET) Update

Posting Date: May 25, 2010

To download the latest update to Version 1.03C, the stand-alone version of the Lab Data Entry Tool (LDET), please click on this link:

<http://www.cdph.ca.gov/programs/aids/Documents/eharsLHJLab.zip>

This LDET version installs to the C: drive of your HARS PC. This update requires a previous installation of LDET. If you have not installed the previous version of LDET (Ver. 1.1) into the encrypted HARS PC, then you must install the previous version first before installing this update.

The name of the file to be downloaded is eharsLHJLAB.ZIP and contains the front-end of LDET called LHJLAB.MDE and a PDF document, Revisions to LDET 1.03C. The front-end of an application consists of the user interface (data entry forms, VBA code, queries, macros, etc.) while the back-end is composed of the database where the case record and lab info are stored. After downloading eharsLHJLAB.ZIP onto your network station's desktop, please copy it to a floppy, or CD, or USB device and transfer it to the desktop of your HARS PC. From the HARS desktop, right-click on eharsLHJLAB.ZIP and extract LHJLAB.MDE to C:\LHJLAB folder in the C: drive of your PC and overwrite the current LHJLAB.MDE file. It will only overwrite the current front-end of the application and will not affect the back-end, or the database of your current LDET.

A new user guide for LDET will be provided in a couple of weeks.

Please contact Victor Borromeo to request for a full installation CD for the new version of LDET if you did not receive the previous version on a CD.

Person to contact for questions or assistance:

Victor Borromeo

Email: victor.borromeo@cdph.ca.gov

Tel. (916) 449-5878