

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Baxter Healthcare Corporation Operating table, Stationary, floor mounting, mobile, and stationary columns, and TruSystem

Recall Date	Product Description	Recalling Firm	Recall Reason
03/21/2025	Baxter Operating table column TS7500 MOBIUS Product Code: 1704695 Baxter Stationary column TruSystem 7500 Product Code: 1717020 Baxter Floor mounting column TruSystem 7500 Product Code: 1717021 Baxter Mobile column TruSystem 7500 Product Code: 1717023 Baxter Mobile column TruSystem 7500 U Product Code: 1730720 Baxter Stationary column TruSystem 7500 U Product Code: 1730731 Baxter Floor mounting column TS 7500 U Product Code: 1730732 Baxter TruSystem 7500 Hybrid (MC) Product Code: 1773204	Baxter Healthcare Corporation	There is a software issue which causes the upper back section to not be operable/adjustable when the "emergency mode" function is enabled.

	<p>Baxter TruSystem 7500 Hybrid (FC) Product Code: 1854085</p> <p>Baxter TruSystem 7500 Hybrid (SC) Product Code: 1854086</p> <p>Baxter TruSystem 7500 Hybrid Plus (FC) Product Code: 1854087</p> <p>Baxter TruSystem 7500 Hybrid Plus (SC) Product Code: 1854088</p> <p>Baxter TruSystem 7500 Hybrid MR IMRIS Product Code: 2067886</p> <p>Baxter TruSystem 7500 Product Code: 4091000</p>		
--	---	--	--

Recall Class	Product Identification	Distribution	Affected Dates
II	<p>UID/DI: 00887761974425 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974418 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974401 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974395 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974388 All serial numbers manufactured until 11/14/2024</p>	11 units in California	February and prior

	<p>UID/DI: 00887761974371 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974364 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974326 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974081 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974074 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974067 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761974050 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761973763 All serial numbers manufactured until 11/14/2024</p> <p>UID/DI: 00887761968639 All serial numbers manufactured until 11/14/2024</p>		
--	--	--	--

For additional information, please visit: [FDA Website \(Operating table column TS7500\)](#), [FDA Website \(Stationary column 7500\)](#), [FDA Website \(Floor mounting column 7500\)](#), [FDA Website \(Mobile column 7500\)](#), [FDA Website \(Mobile column 7500 U\)](#), [FDA Website \(Stationary column 7500 U\)](#), [FDA Website \(Floor mounting column 7500 U\)](#), [FDA Website \(Hybrid \(MC\)\)](#), [FDA Website \(Hybrid \(FC\)\)](#), [FDA Website \(Hybrid \(SC\)\)](#), [FDA Website \(Hybrid Plus \(FC\)\)](#), [FDA Website \(Hybrid Plus \(SC\)\)](#), [FDA Website \(Hybrid MR IMRIS\)](#), and [FDA Website \(TruSystem 7500\)](#).