

## State of California—Health and Human Services Agency

## California Department of Public Health



Food and Drug Branch - Device Recalls

## CALIFORNIA DEVICE RECALL INFORMATION SHEET

Baxter Healthcare Recalls Life2000 Ventilator Series Products for Desaturation Events When Connected to an Oxygen Concentrator

Recall Date	Product Description	Recalling Firm	Recall Reason
3/8/2023	Life2000 Ventilator System  a) REF BT-20-0002, b) REF BT-20-0002A, ventilatory support for the care of individuals	Baxter Healthcare Corporation Deerfield, Illinois	There is the potential for patient desaturation events to occur under certain scenarios when the Life 2000 ventilation system is connected with an oxygen concentrator.
3/8/2023	Life2000 Ventilator PKG, Hospital a) REF BT20007, b) REF BT- 20-0007, ventilatory support for the care of individuals	Baxter Healthcare Corporation Deerfield, Illinois	Potential for patient desaturation events.
3/8/2023	Life2000 Ventilator PKG  a) REF BT200011, b) REF BT- 20-0011, ventilatory support for the care of individuals	Baxter Healthcare Corporation Deerfield, Illinois	Potential for patient desaturation events.
3/8/2023	Life2000 Ventilator  REF MS-01-0118, ventilatory support for the care of individuals	Baxter Healthcare Corporation Deerfield, Illinois	Potential for patient desaturation events.

Recall Class	Product Identification	Distribution	Affected Dates
I	a) REF BT-20-0002, UDI 00887761978089, Serial Numbers: 120370002037 – 1220500011	4279 units Nationwide	January 2023 and Prior

I	a) REF BT20007, UDI 00887761978089, Serial Numbers: L191001-005 - L2005055-898	135 units Nationwide	January 2023 and Prior
I	a) REF BT200011, UDI 00887761978089, Serial Numbers: L190606-014, L191011-040, L200127-015, L200204-013; b) REF BT-20-0011, UDI 00887761978089, Serial Numbers: L200226-012, L200226-013, L200226-015		January 2023 and Prior
I	UDI 00887761978089, Serial Numbers: L1912010-050 - L1912025-007	71 units Nationwide	January 2023 and Prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE <u>FDA WEBSITE</u>

