



California Device Recall Information Sheet

Food and Drug Branch - Device Recalls

DATEX--OHMEDA. INC. Recalls GE Heathcare Giraffe

Recall Date	Product Description	Recalling Firm	Recall Reason
10/17/2024	l '	INC.	GE HealthCare has become aware that for certain Giraffe OmniBed and Giraffe OmniBed Carestation devices, the screw that secures the doors that cover the warmer heaters may not have been torqued to specification. This may lead to the doors becoming loose.

Recall Class	Product Identification	Distribution	Affected Dates
I	All models	Worldwide distribution.	September 2024 and Prior

For additional information, please visit the <u>FDA Website</u>.