

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Draeger Medical Systems, Inc. Recalls Infinity Central Station. Model Number: MS26800

Recall Date	Product Description	Recalling Firm	Recall Reason
1/30/2025	Infinity Central Station. Model Number: MS26800 The Infinity CentralStation (ICS) is intended for centralized monitoring of adult, pediatric and neonatal patient data within the hospital or clinical environment. Centralized monitoring involves the display and management of data from networked patient monitors including the annunciation of visual and audible physiologic parameter alarms at a central monitoring workstation.	Draeger Medical Systems, Inc.	Edifier speakers (shipped in or after April 2023) may not annunciate when configured in analog mode which can lead to delayed treatment to the patient.

Recall Class	Product Identification	Distribution	Affected Dates
II	Infinity CentralStation (ICS) Model Number: MS26800. UDI-DI Number: 04049098001878. Edifier Speaker Model Number: MS34036. UDI-DI Number: 04048675548027. Serial Numbers: 12946M00269, 12946M00314, 12946M00506, 12946M00519, 12946M00518, 12946M00553, 1OQ7WR-00103, 1OQ8CA-00249, 1OQ8CA-00254, 1OQ8CA-00229, 1OQ8CA-00221, 1OQ8CA-00282, 1OQ8CA-00314, 1OQ8CA-00193, 1OQ8CA-00200, 1OQ8CA-00194, 1OQ8CA-00181, 1OQ8CA-00178, 1OQ8CA-00248,	2 devices were distributed to hospitals in California	January 2025 and Prior

10Q8CA-00125, 10Q8CA-00071, 10Q8CA-00032, 10Q8CA-00171, 10Q8CA-00164, 10Q8CA-00149, 10Q8CA-00069, 10Q8CA-00057, 10Q8CA-00153, 10Q8CA-00150, 10Q8CA-00148, 10Q8CA-00165, 10Q8CA-00155, 10Q8CA-00097, 10Q8CA-00113, 10Q8CA-00106, 10Q8CA-00075, 10Q8CA-00052, 10Q8CA-00067, 10Q8CA-00128, 10Q8CA-00132, 10Q8CA-00268, 10Q8CA-00198, 10Q8CA-00255, 10Q8CA-00151, 10Q8CA-00342, 10Q8CA-00391, 10Q8CA-00321, 10Q8CA-00402, 10Q8CA-00395, 10Q8CA-00297, 10Q8CA-00240, 10Q89M-00691, 10Q89M-00690, 10Q89M-00708, 10Q89M-00693, 10Q89M-00687, 10Q89M-00707, 10Q89M-00703, 10Q89M-00679, 10Q89M-00683, 10Q89M-00692, 10Q89M-00610, 10Q89M-00613, 10Q89M-00669, 10Q89M-00635, 10Q89M-00664, 10Q89M-00697, P10Q84F-00080, 10Q89M-00684, 10Q89M-00676, 10Q89M-00695, 10Q8CA-00086, 10Q89M-00677, 10Q8CA-00031, 10Q8CA-00027, 10Q8CA-00065, 10Q8CA-00092, 10Q8CA-00082, 10Q89M-00602, 10Q89M-00680, 10Q89M-00675, 10Q89M-00593, 10Q89M-00686, 10Q89M-00696, 10Q89M-00681, 10Q89M-00704, 10Q89M-00674, 10Q89M-00661, 10Q89M-00663, P10Q89M-00559, P10Q89M-00540, 10Q8CA-00134, 10Q8CA-00095, 10Q8CA-00043, 10Q8CA-00074, 10Q8CA-00762, 10Q8CA-00688, 10Q8CA-00670, 10Q8CA-00637, 10Q8CA-00633, 10Q8E8-00037, 10Q8E8-00053, 10Q8CA-00773, 10Q8CA-00683, 10Q8CA-00500, 10Q8CA-00476, 10Q8CA-00449, 10Q8CA-00420, 10Q8CA-00460, 10Q8CA-00463, 10Q8CA-00507, 10Q8CA-00350, 10Q8CA-00442, 10Q8CA-00433,		
---	--	--

	10Q8CA-00447, 10Q8CA-00601, 10Q8CA-00604, 10Q8CA-00695, 10Q8CA-00513, 10Q8CA-00575, 10Q8CA-00577, 10Q8CA-00591, 10Q8CA-00197, 10Q8CA-00516, 10Q8CA-00541, 10Q8CA-00542, 10Q8CA-00568, 10Q8CA-00499, 10Q8CA-00498, P12051M00227, 10Q8CA-00593, 12946M00532, 12946M00256, 12946M00278, 10Q7WR-00203		
--	---	--	--

For additional information, please visit the [FDA Website](#).