

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Enchroma Inc EnChroma Safety Glasses with Rx Indoor Universal Lenses

Recall Date	Product Description	Recalling Firm	Recall Reason
11/19/2024	<p>EnChroma Safety Glasses with Rx Indoor Universal Lenses:</p> <p>Martinez Cx1 Indoor Rx SKU: Cx.PC.IN.MTZ.BK.Rx</p> <p>Summit Indoor Rx SKU: Cx.PC.IN.SUM.BK.Rx</p> <p>Martinez Black Cx Rx SKU: Cx-3.PC.PT.MTZ.BK.Rx</p>	Enchroma Inc	Safety glasses with Rx indoor universal lenses fail to meet high impact safety specifications, which may not provide adequate protection against impact, posing a risk of eye injury.

Recall Class	Product Identification	Distribution	Affected Dates
II	<p>Order Numbers: #S166556, #S165570, #S166008, #S166145, #S166935, #S167024, #S168061, #S168409, #S168462, #S169061, #S169266, #S170297, #S170428, #S170496, #S171343, #S172249, #S172993, #S173215, #S173512, #S173548, #S173619, #S173729, #S173920, #S174701, #S174818, #S175254, #S175453, #S175681, #S175744, #S176183, #S176257, #S176715, s175356, s163884, #S163421, #s164055, s164308, s164426, s165796, s164324, s162314, s160669, s160666,</p>	2 units in California.	September 6 th and onwards.

	s159817, s161084, #s160752, ecp11546,#ECP11807, #ECP11063, ecp11867, ecp10693, s166343, s174985		
--	--	--	--

For additional information, please visit the [FDA Website](#).