



Food and Drug Branch – Device Recalls

CALIFORNIA DEVICE RECALL INFORMATION SHEET

Invacare recalls TDX SP2 Powered Wheelchair for Wiring Harness Problem

Recall Date	Product Description	Recalling Firm	Recall Reason
3/30/2022	Invacare TDX SP2 Model Nos. TDXSP2 TDXSP2-MCG TDXSP2V	Invacare Corporation Elyria, Ohio	The battery wiring harness on affected wheelchairs may become disconnected during shipment from the factory or during wheelchair use traversing large obstacles or after a jolt or impact to the wheelchair. The disconnection of the battery harness results in loss of power to the wheelchair.

Recall Class	Product Identification	Distribution	Affected Dates
II	UDI (DI): Model Number Device Identifier TDXSP2 00841447100911 TDXSP2-MCG 00841447100966 TDXSP2V 00841447101512 Serial Numbers: 22CE000031 - 22CE001076	6 units in California	March 2022 and Prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch
 MS 7602 • P.O. Box 997435 • Sacramento, CA 95899-7435
 (916) 650-6500 • (916) 650-6650 FAX
 Internet Address: www.cdph.ca.gov

