

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Medline Industries, LP-Northfield Arterial Line Products

Recall Date	Product Description	Recalling Firm	Recall Reason
1/29/2025	Medline General Surgery Tray kits: Arterial Line Tray, Item Number ART600; A line tray, Item Number ART690	Medline Industries, LP - Northfield	The catheter subassembly in the Integrated Arterial Catheter was manufactured with excess material on the catheter hub.
1/29/2025	Medline General Surgery Tray Kits: Arterial line insertion kit, Item Number ART1185B; Arterial line bundle, Item Number ART255; Arterial line tray, Item Number ART350; Arterial line bundle, Item Number ART355; Arterial line bundle, Item Number ART355; Arterial line insertion tray, Item Number ART475; Arterial line bundle - or, Item Number ART535A; Arterial line tray, Item Number ART545A; Arterial line insertion kit, Item Number ART775C; Arterial line tray, Item Number ART840; Arterial line kit, Item Number ART890B	Medline Industries, LP - Northfield	The catheter subassembly in the Integrated Arterial Catheter was manufactured with excess material on the catheter hub.

Recall Class	Product Identification	Distribution	Affected Dates
I	ART600 UDI-DI: 10653160312727 (each), 00653160312720 (case), Lot Number 2024110190; 2) ART690 UDI-DI: 10653160318019 (each), 00653160318012 (case),	58 kits US Nationwide	December 2024 and Prior

	Lot Number 2024092490		
I	ART1185B UDI-DI: 10653160379225(each), 00653160379228(case), Lot Number: 2024102590; ART255 UDI-DI: 10653160367598(each), 00653160367591(case), Lot Number: 2024082690; ART255 UDI-DI: 10653160367598(each), 00653160367591(case), Lot Number: 2024101590; ART350 UDI-DI: 10653160362050(each), 00653160362053(case), Lot Number: 2024093090; ART355 UDI-DI: 10653160361817(each), 00653160361810(case), Lot Number: 2024081290; ART355 UDI-DI: 10653160361817(each), 00653160361810(case), Lot Number: 2024112590; ART475 UDI-DI: 10653160367659(each), 00653160367652(case), Lot Number: 2024100790; ART475 UDI-DI: 10653160367659(each), 00653160367652(case), Lot Number: 2024073090; ART475 UDI-DI: 10653160367659(each), 00653160367652(case), Lot Number: 2024081590; ART475 UDI-DI: 10653160367659(each), 00653160367652(case), Lot Number: 2024090990; ART535A UDI-DI: 10653160367604(each), 00653160367607(case), Lot Number: 2024100890; ART535A UDI-DI: 10653160367604(each), 00653160367607(case), Lot Number: 2024080190; ART545A UDI-DI: 10653160375135(each), 00653160375138(case), Lot Number: 2024102890; ART775C UDI-DI: 10653160366102(each), 00653160366105(case), Lot Number: 2024083090; ART840 UDI-DI: 10653160366966(each), 00653160366969(case), Lot Number: 2024082390; ART840 UDI-DI: 10653160366966(each),	942 kits US Nationwide	December 2024 and Prior

	00653160366969(case), Lot Number: 2024093090; ART890B UDI-DI: 10653160360834(each), 00653160360837(case), Lot Number: 2024100990;		
--	---	--	--

For additional information, please visit this [FDA Website Arterial Line Tray \(ART600\) and A-Line Tray \(ART690\)](#) and this [FDA Website Arterial Line Insertion Kits \(ART1185B, ART775C\), Arterial Line Bundles \(ART255, ART355, ART535A\), Arterial Line Trays \(ART350, ART475, ART545A, ART840, ART890B\)](#).