

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Medline Industries Medline Kits containing Codman Disposable Perforator

Recall Date	Product Description	Recalling Firm	Recall Reason
05/30/2025	Medline Kits containing Codman Disposable Perforator 14 mm: CRANI ADD A PACK SKU: DYNJ30926O CRANIOTOMY SKU: CDS983611N CRANIOTOMY CDS SKU: CDS983467J CRANIOTOMY PACK SKU: DYNJ85807A, DYNJ85927 MAJOR CRANIOTOMY PACK SKU: DYNJ82007A CRANIOTOMY SKU: DYNJ35349C CRANI ACCESSORY PACK SKU: DYNJ59270 CRANI PACK-LF SKU: DYNJ44805M CRANIOTOMY SKU: CDS982719V, CDS982719W, DYNJ904168F	Medline Industries, LP - Northfield	Affected kits contain recalled Integra Lifesciences Codman Disposable Perforators due to inadequate welds on specific lots. Perforators with inadequate welds may break or separate.

	DYNJ908404B, DYNJ908723A, DYNJ908723C CRANIOTOMY #IMF 56646-LF SKU: DYNJVB91001A CRANIOTOMY CDS SKU: CDS981888X, CDS983467I CRANIOTOMY PACK-LF SKU: DYNJ09882G, DYNJ09882I CRANIOTOMY S NASSAU SKU: DYNJ905879G MHC CRANIOTOMY SKU: DYNJ903713K, DYNJ903713L PAD CRANIOTOMY SKU: DYNJ902149J		
--	---	--	--

Recall Class	Product Identification	Distribution	Affected Dates
I	CRANI ADD A PACK SKU: DYNJ30926O UDI-DI: 10195327539153 (ea), 40195327539154 (case) Kit Lots: 24AMF811 CRANIOTOMY SKU: CDS983611N UDI-DI: 10195327466626 (ea), 40195327466627 (case) Kit Lots: 23KBA720 - 23LBJ586 CRANIOTOMY CDS SKU: CDS983467J UDI-DI: 10195327561116 (ea), 40195327561117 (case) Kit Lots: 24ADA818 - 24IDB043	1,853 units nationwide	April and prior

	CRANIOTOMY PACK SKU: DYNJ85807A UDI-DI: 10195327637392 (ea), 40195327637393 (case) Kit Lots: 24CBB843 - 24GBY029 SKU: DYNJ85927 UDI-DI: 10195327575571 (ea), 40195327575572 (case) Kit Lots: 24BMD310, 24EMI836 MAJOR CRANIOTOMY PACK SKU: DYNJ82007A UDI-DI: 10195327635466 (ea), 40195327635467 (case) Kit Lots: 24CMA496 - 24HMI978 CRANIOTOMY SKU: DYNJ35349C UDI-DI: 10888277858060 (ea), 40888277858061 (case) Kit Lots: 24CBJ925 - 24FBO318 CRANI ACCESSORY PACK SKU: DYNJ59270 UDI-DI: 10888277771512 (ea), 40888277771513 (case) Kit Lots: 23JDB683 - 24HDA822 CRANI PACK-LF SKU: DYNJ44805M UDI-DI: 10195327506872 (ea), 40195327506873 (case) Kit Lots: 23JBI747 - 24GBJ699 CRANIOTOMY SKU: CDS982719V UDI-DI: 10195327510831 (ea), 40195327510832 (case) Kit Lots: 23KBP461 - 24BBI145 SKU: CDS982719W UDI-DI: 10195327645090 (ea), 40195327645091 (case) Kit Lots: 24DBS813 SKU: DYNJ904168F UDI-DI: 10195327395728 (ea), 40195327395729 (case) Kit Lots: 24DBU188		
--	---	--	--

SKU: DYNJ908404B UDI-DI: 10195327410506 (ea), 40195327410507 (case) Kit Lots: 23IDB680 - 24IDA840 SKU: DYNJ908723A UDI-DI: 10195327403270 (ea), 40195327403271 (case) Kit Lots: 23IBO314 - 24ABK090 SKU: DYNJ908723C UDI-DI: 10198459009785 (ea), 40198459009786 (case) Kit Lots: 24IBM205 - 25ABM356 CRANIOTOMY #IMF 56646-LF SKU: DYNJVB91001A UDI-DI: 10195327466466 (ea), 40195327466467 (case) Kit Lots: 23JLA765 - 24ELA325 CRANIOTOMY CDS SKU: CDS981888X UDI-DI: 10195327144760 (ea), 40195327144761 (case) Kit Lots: 24ABT432 - 25CBF098 SKU: CDS983467I UDI-DI: 10195327457945 (ea), 40195327457946 (case) Kit Lots: 23KDA336 CRANIOTOMY PACK-LF SKU: DYNJ09882G UDI-DI: 10195327574055 (ea), 40195327574056 (case) Kit Lots: 24CBF778 SKU: DYNJ09882I UDI-DI: 10198459058172 (ea), 40198459058173 (case) Kit Lots: 24HBN426 CRANIOTOMY S NASSAU SKU: DYNJ905879G UDI-DI: 10195327512996 (ea), 40195327512997 (case) Kit Lots: 23KBI086 - 24FBE409 MHC CRANIOTOMY SKU: DYNJ903713K UDI-DI: 10195327201579 (ea),		
---	--	--

	40195327201570 (case) Kit Lots: 23JBK426 SKU: DYNJ903713L UDI-DI: 10195327531218 (ea), 40195327531219 (case) Kit Lots: 23KBT939 - 24CBH771 PAD CRANIOTOMY SKU: DYNJ902149J UDI-DI: 10195327554323 (ea), 40195327554324 (case) Kit Lots: 23KBP275 - 24IBN443		
--	---	--	--

For additional information, please visit: [FDA Website \(1\)](#), [FDA Website \(2\)](#), and [FDA Website \(3\)](#).