

California Device Recall Information Sheet

Food and Drug Branch – Device Recalls

Medtronic MiniMed, Inc. MiniMed and Paradigm Insulin Pump

Recall Date	Product Description	Recalling Firm	Recall Reason
03/03/2025	MiniMed insulin pump REF: 620G (MMT-1750) MiniMed insulin pump REF: 630G (MMT-1714, MMT-1715, MMT-1754, MMT-1755) MiniMed insulin pump REF: 640G (MMT-1711, MMT-1712, MMT-1751, MMT-1752) MiniMed insulin pump REF: 670G (MMT-1760, MMT-1761, MMT-1762, MMT-1780, MMT-1781, MMT-1782) MiniMed insulin pump REF: 700G (MMT-1801, MMT-1805, MMT-1850, MMT-1851) MiniMed insulin pump REF: 720G (MMT-1809, MMT-1810, MMT-1859, MMT-1860, MMT-1867) MiniMed insulin pump REF: 740G (MMT-1811, MMT-1812, MMT-1861, MMT-1862) MiniMed insulin pump REF: 770G (MMT-1880, MMT-1881, MMT-1882, MMT-1890,	Medtronic MiniMed, Inc.	All Medtronic MiniMed pumps (Paradigm series, Minimed 530G, 6XX series, and 7XX series) were found to have the potential for abnormal insulin delivery due to changes in air pressure, which can be found during airplane take-off or landing. An unexpected additional dose of insulin, which can occur during take-off, can result in severe hypoglycemia, altered mental status, seizure, coma, or death. An unexpected under-delivery of insulin, which can occur during landing, can result in severe hyperglycemia, dehydration, diabetic ketoacidosis, or death.

	MMT-1891, MMT-1892) MiniMed insulin pump REF: 780G (MMT-1884, MMT-1885, MMT-1886, MMT-1894, MMT-1895, MMT-1896) Paradigm insulin pump REF: MMT-522, MMT-523, MMT-551, MMT-722, MMT-723, MMT-751 Paradigm REF: MMT-712 and MMT-715 Paradigm insulin pump REF: MMT-754		
--	---	--	--

Recall Class	Product Identification	Distribution	Affected Dates
II	REF: 620G (MMT-1750) All Lot Numbers All UDI numbers REF: 630G (MMT-1714, MMT-1715, MMT-1754, MMT-1755) All Lot Numbers All UDI numbers REF: 640G (MMT-1711, MMT-1712, MMT-1751, MMT-1752) All Lot Numbers All UDI numbers REF: 670G (MMT-1760, MMT-1761, MMT-1762, MMT-1780, MMT-1781, MMT-1782) All Lot Numbers All UDI numbers REF: 700G (MMT-1801, MMT-1805, MMT-1850, MMT-1851) All Lot Numbers All UDI numbers	Nationwide	January and prior

	REF: 720G (MMT-1809, MMT-1810, MMT-1859, MMT-1860, MMT-1867) All Lot Numbers All UDI numbers REF: 740G (MMT-1811, MMT-1812, MMT-1861, MMT-1862) All Lot Numbers All UDI numbers REF: 770G (MMT-1880, MMT-1881, MMT-1882, MMT-1890, MMT-1891, MMT-1892) All Lot Numbers All UDI numbers REF: 780G (MMT-1884, MMT-1885, MMT-1886, MMT-1894, MMT-1895, MMT-1896) All Lot Numbers All UDI numbers REF: MMT-522, MMT-523, MMT-551, MMT-722, MMT-723, MMT-751 All Lot Numbers All UDI numbers REF: MMT-712 and MMT-715 All Lot Numbers All UDI numbers REF: MMT-754 All Lot Numbers All UDI numbers		
--	--	--	--

For additional information, please visit: [FDA Website \(620G\)](#), [FDA Website \(630G\)](#), [FDA Website \(640G\)](#), [FDA Website \(670G\)](#), [FDA Website \(700G\)](#), [FDA Website \(720G\)](#), [FDA Website \(740G\)](#), [FDA Website \(770G\)](#), [FDA Website \(780G\)](#), [FDA Website \(Paradiqm MMT-522\)](#), [FDA Website \(Paradiqm MMT-712\)](#), and [FDA Website \(Paradiqm MMT-754\)](#).