



# California Medical Device Recall Information



## Recall Name

### Maquet Medical Systems Recalls Flow-i Anesthesia Systems Due to a Loose Patient Cassettes

Recall Date	Product Description	Recalling Firm	Recall Reason
04/02/15	Flow-i Anesthesia Systems	<b>Maquet Medical Systems</b> Wayne, NJ	<i>Patient cassettes, which are the center of gas flow in the system, have come loose.</i>  <i>This may cause anesthesia gas to leak and could prevent the ventilator from providing breathing support if not corrected immediately.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	<b>Article Numbers:</b>  <ul style="list-style-type: none"> <li>• 66 77 200</li> <li>• 66 77 300</li> <li>• 66 77 400</li> </ul>	CA, nationwide	Manufacturing and Distribution Dates:  July 7, 2011 to May 5, 2014

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm453492.htm>