



California Medical Device Recall Information



Recall Name

**Maquet Medical Systems Recalls TigerPaw System II Surgical Staple –
May Cause Tears and Bleeding in Heart Tissue**

Recall Date	Product Description	Recalling Firm	Recall Reason
03/25/15	TigerPaw System II	Maquet Medical Systems Wayne, NJ	<i>Potential for incomplete closure of the surgical staple may result in tissue tears and/or bleeding.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	All Lots Recalled Part Number: <ul style="list-style-type: none"> • C-TP-1507 • C-TP-1509 	CA , nationwide	Distributed: April 1, 2013 through March 23, 2015

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm446087.htm>