



California Medical Device Recall Information



Recall Name

Verathon Inc. Recalls GlideScope Laryngoscope Blades Due to Breaking Blades

Recall Date	Product Description	Recalling Firm	Recall Reason
10/25/12	GlideScope GVL Video Laryngoscope Reusable Blades	Verathon, Inc. Bothell, WA	<i>Prone to developing cracks and/or breaking across the tip of the blade</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	<ul style="list-style-type: none"> GlideScope GVL 3 Model # 0574-0007 S/Ns: MD10500 to MD112387 GlideScope GVL 4 Model #0574-0001 S/Ns: LG105000 to LG112758 GlideScope GVL 5 Model #0574-0030 S/Ns: XL105000 to XL111798 <p>The reusable blades can be identified by the serial number engraved on the label found on the handle of the blade.</p>	CA, nationwide	<p>Manufactured from December 1, 2010 to August 31, 2011;</p> <p>Distributed from December 14, 2010 to March 6, 2012</p>

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/Safety/Recalls/ucm330305.htm>