



California Medical Device Recall Information



Recall Name

Western/Scott Fetzer Company Recalls OxyTote Portable Oxygen Unit Due to Potential for the Cylinder to Ignite

Recall Date	Product Description	Recalling Firm	Recall Reason
1/31/15	OxyTote Portable Oxygen Unit	Western/Scott Fetzer Company Westlake, OH	<i>The company received reports that when the OxyTote is mishandled or dropped, the oxygen cylinder may ignite causing an internal flash fire and the canister to burst.</i>
Recall Class	Product Identification	Distribution	Affected Dates
I	Affected Lot and Model Numbers	CA, nationwide	Manufacturing and Distribution Dates: January 1, 2009 to September 30, 2014

FOR ADDITIONAL INFORMATION, PLEASE VISIT:

<http://www.fda.gov/MedicalDevices/Safety/ListofRecalls/ucm458740.htm>