



**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**Baxter recalls Inline Ventilator Adaptor for Risk of Reduced Ventilation or Pneumothorax**

Recall Date	Product Description	Recalling Firm	Recall Reason
4/26/2022	<b>In-Line ventilator adaptor</b>	<b>Baxter Healthcare Corporation</b> Deerfield, Illinois	There is a potential of reduced oxygenation or pneumothorax /barotrauma to occur when the Volara system is used with in-line ventilator in home care environment.
4/26/2022	<b>In-Line ventilator adaptor</b>	<b>Baxter Healthcare Corporation</b> Deerfield, Illinois	Potential of reduced oxygenation or pneumothorax.

Recall Class	Product Identification	Distribution	Affected Dates
I	<b>M08473 OPTIMUS OLE AC PAT.CIRCUIT KIT</b> UDI-DI: 00887761981492.	259 Units Nationwide including California	February 2020 to April 2022
I	<b>M07937 MODULE, OPTIMUS HANDSET 2;</b> UDI-DI: 00887761984622.	9 Units Nationwide including California	October 2019 to April 2022

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

