



## CALIFORNIA DEVICE RECALL INFORMATION SHEET

**Caire, Inc. Recalls CAIRE Liberator Liquid Oxygen System Units  
due to inconsistencies in weld penetration on the inner bottle**

Recall Date	Product Description	Recalling Firm	Recall Reason
12/14/2022	<b>CAIRE Liberator 45, MODEL 13261699, Liquid Oxygen System Unit</b> 10 devices	<b>Caire, Inc.,</b> Ball Ground, Georgia	An audit discovered some inconsistencies in weld penetration on the longitudinal weld seam on the inner bottle of the device on some of the units.
12/14/2022	<b>CAIRE Liberator 45, MODEL 13262253, Liquid Oxygen System Unit</b> 70 Devices	<b>Caire, Inc.,</b> Ball Ground, Georgia	An audit discovered some inconsistencies in weld penetration on the longitudinal weld seam on the inner bottle of the device on some of the units.
12/14/2022	<b>CAIRE Liberator 30, MODEL 13337403, Liquid Oxygen System Unit</b> 5 Devices	<b>Caire, Inc.,</b> Ball Ground, Georgia	An audit discovered some inconsistencies in weld penetration on the longitudinal weld seam on the inner bottle of the device on some of the units.
12/14/2022	<b>CAIRE Liberator 20, MODEL 13256195, Liquid Oxygen System Unit</b> 3 Devices	<b>Caire, Inc.,</b> Ball Ground, Georgia	An audit discovered some inconsistencies in weld penetration on the longitudinal weld seam on the inner bottle of the device on some of the units.

Recall Class	Product Identification	Distribution	Affected Dates
I	UDI/DI M766132616990, Serial Numbers: CBB3022310887, CBB3022310888, CBB3022310852, CBB3022310823, CBB3022310792, CBB3022310843, CBB3022310844, CBB3022310827, CBB3022310842, CBB3022310845	Worldwide Distribution including California	October 2022 and prior
I	UDI/DI M766132622530, Serial Numbers: CBB3022311044, CBB3022311042, CBB3022311043, CBB3022311040, CBB3022311242, CBB3022311098, CBB3022311093, CBB3022311094, CBB3022311077, CBB3022311053, CBB3022311049, CBB3022311065, CBB3022311051, CBB3022311075, CBB3022311078, CBB3022311097, CBB3022311096, CBB3022311079, CBB3022311045, CBB3022311046, CBB3022311236, CBB3022311230, CBB3022311217, CBB3022311231, CBB3022311164, CBB3022311225, CBB3022311226, CBB3022311133, CBB3022311192, CBB3022311163, CBB3022311216, CBB3022311228, CBB3022311227, CBB3022311162, CBB3022311198, CBB3022311211, CBB3022311202, CBB3022311155,	Worldwide Distribution including California	October 2022 and prior

	CBB3022311135, CBB3022311171, CBB3022311121, CBB3022311106, CBB3022311109, CBB3022311108, CBB3022311084, CBB3022311083, CBB3022311092, CBB3022311081, CBB3022311082, CBB3022311141, CBB3022311131, CBB3022311136, CBB3022311183, CBB3022311178, CBB3022311130, CBB3022311179, CBB3022311139, CBB3022311177, CBB3022311105, CBB3022311146, CBB3022311172, CBB3022311140, CBB3022311125, CBB3022311169, CBB3022311068, CBB3022311176, CBB3022311170, CBB3022311180, CBB3022311057, CBB3022311129		
I	UDI/DI M766133374030, Serial Numbers: CBB3022300368, CBB3022300390, CBB3022300470, CBB3022300472, CBB3022300479	Worldwide Distribution including California	October 2022 and prior
I	UDI/DI M766133374030M7661325619 50, Serial Numbers: CBB3022320073, CBB3022320077, CBB3022320074	Worldwide Distribution including California	October 2022 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch  
MS 7602 ● P.O. Box 997435 ● Sacramento, CA 95899-7435  
(916) 650-6500 ● (916) 650-6650 FAX  
Internet Address: [www.cdph.ca.gov](http://www.cdph.ca.gov)

