



Food and Drug Branch – Device Recalls

**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**Draeger Medical Recalls Ventstar Coax 180 For Loosening Of Glued Connections**

Recall Date	Product Description	Recalling Firm	Recall Reason
5/24/2023	<b>Ventstar Coax 180</b> MP03374	<b>DRAEGER MEDICAL</b> Telford, Pennsylvania	Glued connections of the breathing circuit can become loose before or during ventilation resulting in partial or complete detachment of pieces.

Recall Class	Product Identification	Distribution	Affected Dates
I	<b>Ventstar Coax 180</b> UDI-DI: 04048675545576; All lots.	20 Units in California	April 2023 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch  
 MS 7602 • P.O. Box 997435 • Sacramento, CA 95899-7435  
 (916) 650-6500 • (916) 650-6650 FAX  
 Internet Address: [www.cdph.ca.gov](http://www.cdph.ca.gov)

