



**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**NeuroLogica Corporation Recalls Samsung Digital Diagnostic Mobile X-ray System due to the Welding of the Moving Arm Frame**

Recall Date	Product Description	Recalling Firm	Recall Reason
2/21/2024	<b>Samsung Digital Diagnostic Mobile X-ray System</b> Model GM85.	<b>NeuroLogica Corporation</b> Danvers, Massachusetts	This product is being recalled due to an issue with the welding of the moving arm frame. If the arm falls, this poses a risk of bodily harm.

Recall Class	Product Identification	Distribution	Affected Dates
II	<b>Samsung Digital Diagnostic Mobile X-ray System</b> UDI-DI: 08806088582412	73 Units in California California	February, 2024 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

