



## CALIFORNIA DEVICE RECALL INFORMATION SHEET

### Olympus Corporation of the Americas Bronchofiberscope and Bronchovideo Scope for Endobronchial Combustion

Recall Date	Product Description	Recalling Firm	Recall Reason
7/26/2023	<b>Olympus OES BRONCHOFIBERSCOPE BF Type</b> Model Numbers BF-XT40, BF-P60, BF-MP60, BF-1T60	<b>Olympus Corporation of the Americas</b> Center Valley, Pennsylvania	There have been complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.
7/26/2023	<b>Olympus BRONCHOFIBERSCOPE BF Type</b> Model Numbers BF-PE2, BF-TE2	<b>Olympus Corporation of the Americas</b> Center Valley, Pennsylvania	There have been complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.
7/26/2023	<b>Olympus BRONCHOVIDEOSCOPE BF Type</b> Model Numbers BF-P150, BF-1T150	<b>Olympus Corporation of the Americas</b> Center Valley, Pennsylvania	There have been complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.
7/26/2023	<b>Olympus EVIS EXERA BRONCHOVIDEOSCOPE BF TYPE</b> Model Number BF-XT160	<b>Olympus Corporation of the Americas</b> Center Valley, Pennsylvania	There have been complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.
7/26/2023	<b>Olympus EVIS EXERA II</b>	<b>Olympus</b>	There have been

	<b>BRONCHOVIDEOSCOPE BF TYPE</b> Model Numbers BF-P180, BF-Q180, BF-Q180-AC, BF-1T180, BF-1TQ180	<b>Corporation of the Americas</b> Center Valley, Pennsylvania	complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.
7/26/2023	<b>Olympus EVIS EXERA III BRONCHOVIDEOSCOPE BF TYPE</b> Model Numbers BF-H190, BF-Q190, BF-XT190, BF-1TH190	<b>Olympus Corporation of the Americas</b> Center Valley, Pennsylvania	There have been complaints of endobronchial combustion during therapeutic laser procedures with the bronchoscope.

Recall Class	Product Identification	Distribution	Affected Dates
I	<b>Olympus OES BRONCHOFIBERSCOPE BF Type</b> All serial numbers BF-XT40** OES BRONCHOFIBERSCOPE OLYMPUS BF TYPE XT40 UDI: 04953170055980 BF-P60 OES BRONCHOFIBERSCOPE OLYMPUS BF TYPE P60 UDI: 04953170339288 BF-MP60 OES BRONCHOFIBERSCOPE OLYMPUS BF TYPE MP60 UDI: 04953170338394 BF-1T60 OES BRONCHOFIBERSCOPE OLYMPUS BF TYPE 1T60 UDI: 04953170339264	1535 Units Nationwide	June 2023 and Prior
I	<b>Olympus BRONCHOFIBERSCOPE BF Type</b> All serial numbers BF-PE2 BRONCHOFIBERSCOPE OLYMPUS BF TYPE PE2 UDI: 04953170339974 BF-TE2 BRONCHOFIBERSCOPE OLYMPUS BF TYPE TE2 UDI: 04953170339998	4 Units Nationwide	June 2023 and Prior
I	<b>Olympus BRONCHOVIDEOSCOPE BF Type</b>	9 Units Nationwide	June 2023 and Prior

	<p>All serial numbers BF-P150*  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE P150 UDI:  04953170288876 BF-1T150  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE 1T150  UDI: 04953170288968</p>		
I	<p><b>Olympus EVIS EXERA  BRONCHOVIDEOSCOPE BF  TYPE</b>  All serial numbers BF-XT160*  EVIS EXERA  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE XT160  UDI: 04953170340147</p>	284 Units Nationwide	June 2023 and Prior
I	<p><b>Olympus EVIS EXERA II  BRONCHOVIDEOSCOPE BF  TYPE</b>  All serial numbers BF-P180*  EVIS EXERA II  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE P180 UDI:  04953170339288 BF-Q180**  EVIS EXERA II  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE Q180 UDI:  04953170339301 BF-Q180-AC*  EVIS EXERA II  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE Q180-AC  UDI: 04953170340086 BF-  1T180* EVIS EXERA II  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE 1T180  UDI: 04953170339325 BF-  1TQ180* EVIS EXERA II  BRONCHOVIDEOSCOPE  OLYMPUS BF TYPE 1TQ180  UDI: 04953170339349</p>	17,609 Units Nationwide	June 2023 and Prior
I	<p><b>Olympus EVIS EXERA III  BRONCHOVIDEOSCOPE BF  TYPE</b>  All serial numbers BF-H190  EVIS EXERA III  BRONCHOVIDEOSCOPE  OLYMPUS BF-H190 UDI:  04953170434754 BF-Q190  EVIS EXERA III  BRONCHOVIDEOSCOPE</p>	9,684 Units Nationwide	June 2023 and Prior

	OLYMPUS BF-Q190 UDI: 04953170335198 BF-XT190 EVIS EXERA III BRONCHOVIDEOSCOPE OLYMPUS BF-XT190UDI: 04953170402470 BF-1TH190 EVIS EXERA III BRONCHOVIDEOSCOPE OLYMPUS BF-1TH190 UDI: 04953170434778		
--	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--	--

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

---

CDPH Food and Drug Branch  
MS 7602 ● P.O. Box 997435 ● Sacramento, CA 95899-7435  
(916) 650-6500 ● (916) 650-6650 FAX  
Internet Address: [www.cdph.ca.gov](http://www.cdph.ca.gov)

