



Food and Drug Branch – Device Recalls

CALIFORNIA DEVICE RECALL INFORMATION SHEET

Phillips Medical Systems Recalls CombiDiagnost due to potential for the RF to incorrectly display a previous patient’s radiography images.

Recall Date	Product Description	Recalling Firm	Recall Reason
2/21/2024	CombiDiagnost R90 1.1 (709031)	PHILIPS MEDICAL SYSTEMS Cambridge, Massachusetts	This product is being recalled due to there being a potential that the Radio Fluoroscopy (RF) viewer will also display a previous patient's radiography images. If the issue occurs, there will be differences in image content, image format and image size.

Recall Class	Product Identification	Distribution	Affected Dates
II	CombiDiagnost Model: R90 1.1 (709031); UDI DI: 00884838101456; Software version 1.1.1 or higher; Serial Numbers: 10001004, 10001005, 10001006, 10001008, 10001009, 10001010, 10001011, 10001012, 10001013, 10001018, 10001020, 10001022, 10001024, 10001026, 10001027, 10001028, 10001029, 10001030, 10001031, 10001032, 10001035, 10001036, 10001043, 10001044, 10001045, 10001046, 10001047, 10001048, 10001049, 10001050, 10001051, 10001052,	2 Units in California California	February, 2024 and prior

	10001053, 10001054, 10001055, 10001056, 10001057, 10001058, 10001059, 10001060, 10001061, 10001062, 10001063, 10001064, 10001065, 10001066, 10001067, 10001068, 10001069, 10001070, 10001071, 10001072, 10001073, 10001074, 10001075, 10001076, 10001077, 10001078, 10001079, 10001080, 10001081, 10001083, 10001084, 10001085, 10001086, 10001087, 10001088, 10001089, 10001090, 10001091, 10001092, 10001093, 10001094, 10001095, 10001096, 10001097, 10001098, 10001099, 10001100, 10001101, 10001102, 10001103, 10001104, 10001105, 10001106, 10001107, 10001108, 10001109, 10001110, 10001111, 10001112, 10001113, 10001114, 10001115, 10001116, 10001117, 10001118, 10001119, 10001120, 10001121, 10001122, 10001123, 10001124, 10001125, 10001126, 10001127, 10001128, 10001129, 10001130, 10001131, 10001132, 10001133, 10001134, 10001135, 10001136, 10001137, 10001138, 10001139, 10001140, 10001141, 10001142, 10001143, 10001144, 10001145, 10001147, 10001148, 10001149, 10001150, 10001152, 10001201;		
--	---	--	--

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch
MS 7602 • P.O. Box 997435 • Sacramento, CA 95899-7435
(916) 650-6500 • (916) 650-6650 FAX
Internet Address: www.cdph.ca.gov

