



CALIFORNIA DEVICE RECALL INFORMATION SHEET

Resmed Recalls Airtouch And Airfit Masks For Magnetic Interference

Recall Date	Product Description	Recalling Firm	Recall Reason
12/27/2023	Airtouch F20 Full Face Mask and User Guide	RESMED LTD. Bella Vista, Australia	Masks non-invasive interface for channeling airflow to patients have magnets that have potential magnetic interference with active medical implants that interact with magnets, and metallic implants /objects containing ferromagnetic material so contraindications and warnings about safe distance to medical devices /implants are being updated for patients using masks or anyone in close physical contact
12/27/2023	Airtouch N20 Nasal Mask and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference
12/27/2023	Airfit F30 Full Face Mask and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference
12/27/2023	Airfit N10 Nasal Masks and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference
12/27/2023	Airfit F20 Full Face Mask and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference

12/27/2023	Airfit F30i Full Face Mask and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference
12/27/2023	Airfit N20 Nasal Mask and User Guide	RESMED LTD. Bella Vista, Australia	Magnetic Interference

Recall Class	Product Identification	Distribution	Affected Dates
I	Airtouch F20 All mask lots used with User Guide 638220/version 2020-02, and prior. UDI-DI/GTIN: 619498630028 - 619498630370	3,816,138 Units Nationwide	February 2020 and prior
I	Airtouch N20 All mask lots used with User Guide 638218/version 2020-02, and prior. UDI-DI/GTIN: 619498639069 - 619498639212	424309 Units Nationwide	February 2020 and prior
I	Airfit F30 All mask lots used with User Guide 638234/version 2020-07, and prior. UDI-DI/GTIN: 619498641048 - 619498641147	2,660,776 Units Nationwide	July 2020 and prior
I	Airfit N10 All mask lots used with User Guide 628153/version 2020-07, and prior. UDI-DI/GTIN: 619498632008 - 619498632176	199796 Units Nationwide	July 2020 and prior
I	Airfit F20 All mask lots used with User Guide 638220/version 2020-02, 638238/version 2020-12, 638221/version 2020-01, and prior. UDI-DI/GTIN: 619498634019 - 619498634194	11,885,039 Units Nationwide	December 2020 and prior
I	Airfit F30i All mask lots used with User Guide 638225/version 2020-02, and prior. UDI-DI/GTIN:	2037323 Units Nationwide	February 2020 and prior

	619498633326 - 619498633074		
I	Airfit N20 All mask lots used with User Guide 638218/version 2020-02, 638243/version 2020-06, and prior. UDI-DI/GTIN: 619498635016 - 619498635153	7,724,968 Units Nationwide	June 2020 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE FDA WEBSITE at Airtouch[F20](#), Airtouch[N20](#), Airtouch[F30](#), Airfit[N10](#), Airfit[F20](#), Airfit[F30i](#) and Airfit[N20](#)

CDPH Food and Drug Branch
MS 7602 ● P.O. Box 997435 ● Sacramento, CA 95899-7435
(916) 650-6500 ● (916) 650-6650 FAX
Internet Address: www.cdph.ca.gov

