

Micro Learning Video: CCC Specimen Events

In this video, you will learn how to enter or undo a specimen event.

1. Click **Specimens**.
2. Click the required specimen record.
3. From the Summary tab, press **Enter** or click the **Magnifying Glass** icon in the Specimen Event field to view a list of specimen events.
4. Select the required event for that specimen type, either MSAFP or cfDNA.
5. Enter any **Coordinator Notes**, if applicable.
6. Click **Save & Close** or the **Save** button to update the changes.
7. To undo a previous specimen event, select specimen event 99 from the lookup field.

There must be a specimen event saved within the last 3 days to be able to save specimen event 99 or an error message will appear.